



FORSURE

FORSURE Salmonella Antigen Kit are used for the qualitative and semi-quantitative detection of salmonella bacteria agglutinins in bacterial infection. Antibodies are formed in human infection cases with microbiological agents. Mixing these antibodies with the corresponding homologous antigens causes agglutination. The agglutination is macroscopically visible. For screening purposes, the antigens can be used in the qualitative rapid slide test. Positive results can be confirmed with the quantitative tube test to verify the antibody titer.



Cat Nos.	FS-S8Ag-FA5-010
Description	Forsure Salmonella Antigen Kits
Content	5ml OA Reagent, 5ml HA Reagent 5ml OB Reagent, 5ml HB Reagent 5ml OC Reagent, 5ml HC Reagent 5ml OD Reagent, 5ml HD Reagent 1ml Positive Control 1ml Negative Control

TS Diagnostics Sdn Bhd

No 7 & 9, Jalan TPP 5/11,
Taman Perindustrian Puchong,
47160 Puchong, Selangor, Malaysia.

Ph number +6 03 8068 7311
Email contactus@forsure-test.com
Website www.forsure-test.com

